



May 2025

Dear Valued Client,

The Abbott Alinity s and Grifols Panther package inserts for infectious disease markers specify that cadaveric specimens must be collected within 24 hours of the individual's time of death.

To comply with these requirements, the time of death must be clearly noted on the test requisition form when submitting cadaveric samples for infectious disease testing. Any specimens collected more than 24 hours after death will be rejected and testing will not be performed.

We appreciate the opportunity to serve your testing needs. Please contact ViroMed Account Management at 800-582-0077 or ViroMed_Acctmgmt@labcorp.com if you have any questions.

Sincerely,

A handwritten signature in black ink, appearing to read "Ajay Grover", with a long horizontal flourish extending to the right.

Ajay Grover, MS, PhD, BCMAS, HCLD (ABB)
Laboratory Director
ViroMed Laboratories